


# Instructions for use

## **ROOKE® BK Rigid Protector** *with Soft Interface*



1. Open BK Rigid Protector. Center patella in the middle of device.



2. For best protection from falls use two spacers and extend the distal pad if necessary.



3. Place tongue over spacers and residual limb.



4. Trim to fit if necessary. Cut in between dark trim lines.



5. Secure flaps.



6. Buckle distal end. Trim off excess strap if needed.



7. Trim proximal end if necessary. Cut in between the dark trim lines. Add extender if needed.


Continued on back....

## Instructions for use

### **ROOKE® BK Rigid Protector** with Soft Interface



8. Secure flaps.



9. Buckle both thigh buckles. Trim off excess strap if needed.




10. Apply patella pad securely using both hands.



11. Fasten the patella pad.



12. Adjust waist belt & fit to size by trimming off any excess length. Attach hook and loop end to BK protector.



#### Accessory Pack Includes:

- A. Waist Belt.
- B. 2 spacers for residual limb.
- C. Scissors for trimming if necessary.
- D. Extender for either proximal/distal flaps.

 **Osborn Medical**

P.O. Box 324, 100 West Main St.  
Utica, MN 55979

1.800.535.5865 507.932.5028 Fax 507.932.5044

[www.osbornmedical.com](http://www.osbornmedical.com)