





ROOKE® Soft AFO™

Rooke® Soft AFO™ helps maintain Dorsiflexion with the malleable splint providing full plantar support for the affected limb.



The Rooke® Soft AFO™ is constructed by inserting the foam covered aluminum splint into the specially fabricated pocket of the Rooke HFS Boot.*

Once inserted, all hook and loop fasteners are connected and the Rooke® AFO™ is ready for use.



**Please note. Reversing this action allows for the use of the Rooke HFS Boot independent of the AFO.*

Ordering Information

Order No.	Product Description	Size
HFS104	Rooke Soft AFO (includes Rooke HFS Boot with removable padded aluminum foot drop splint)	One size fits all.
Prices subject to change. F.O.B. Utica, MN, U.S.A.		

Osborn Medical

P.O. Box 324, 100 West Main St., Utica, MN 55979
1.800.535.5865 507.932.5028 Fax 507.932.5044
www.RookeBoot.com