

Rooke® Soft AFO™ helps maintain Dorsiflexion with the malleable splint providing full plantar support for the affected limb.



The Rooke® Soft AFO™ is constructed by inserting the foam covered aluminum splint into the specially fabricated pocket of the Rooke HFS Boot.\*

Once inserted, all hook and loop fasteners are connected and the Rooke® AFO™ is ready for use.



*\*Please note. Reversing this action allows for the use of the Rooke HFS Boot independent of the AFO.*

### Ordering Information

Order No.	Product Description	Size
<b>HFS104</b>	Rooke Soft AFO (includes Rooke HFS Boot with removable padded aluminum foot drop splint)	One size fits all.
Prices subject to change. F.O.B. Utica, MN, U.S.A.		

### Osborn Medical

P.O. Box 324, 100 West Main St., Utica, MN 55979  
 1.800.535.5865 507.932.5028 Fax 507.932.5044  
[www.RookeBoot.com](http://www.RookeBoot.com)