

ROOKE® PERIOPERATIVE BOOT

Optimal Insulation and Protection for the Healing Limb



The Rooke® Perioperative Boot

is ideal for patients recovering from surgery, being treated for ischemia, neuropathy or other related conditions. Designed to completely surround and the foot, heel and calf, it is also useful for patients suffering from pressure ulcers.

Gentle warming of the limb

is achieved by reducing heat loss frequently associated with ischemic conditions. The boot does not create heat, but acts as an insulating barrier that retains body-generated heat. The Rooke® Perioperative Boot is super-insulated, containing a 140 gram layer of Thinsulate® plus 1/8th inch of foam padding.

Protection of an ischemic or healing limb is generally regarded as an important aspect of medical management.

Placing such a limb in a Rooke® Perioperative Boot will give it the extra protection it requires.



PERIOPERATIVE BOOT DESIGN FEATURES

Bi-Layer construction

The Rooke® Perioperative Boot has an inner layer of seamless soft fleece. No seams mean that there are no rough edges for sensitive limbs to rub against. The outer layer is constructed from low-friction fabric, which minimizes friction between the boot and bedding. The boot was also designed to be soft-to-the-touch by the patient's other limb.

Molded Foam Cradle

The boot's foam cradle surrounds the dorsal region of the limb, absorbs pressure from bedding, which touches the limb and acts as a barrier against additional trauma. The foam cradle provides support to the heel and leg, while also providing extra protection. The cradle's anatomically correct design distributes weight evenly to prevent pressure sores during prolonged bed rest. The cradle also suspends the limb and limited off-loading is achieved.

Hook and Loop Closures

Adjustable fasteners allow for easy operation by nursing staff or the patient. The four fasteners are strategically placed so the boot may be worn in a variety of positions. This adaptability allows the boot to be manipulated for use with many conditions.



1.800.535.5865
or 507.932.5028
Fax 507.932.5044

100 West Main St., P.O. Box 324
Utica, MN 55979

www.osbornmedical.com